



Gentle 100% Ultrasound Technology for better Oral Health

Emmi-dent Oral Waver (with added gentle vibration feature) Quick Start Guide

Welcome to Your completely new oral care cleaning experience – with 100% Ultrasound waves. Importantly – this is not a regular toothbrush. Cleaning does NOT happen via abrasive brushing but via gentle 100% Ultrasound wave cleaning technology. Please follow usage instructions carefully.

1) **IMPORTANT: Before 1st usage, charge device for at least 24 hours.**

- Place device in charging unit (blue charging LED will turn on) and charge for at least 24 hrs.
- When fully charged and turned on, three green LED lights will light up sequentially.

2) To clean teeth and gums

- Wet the bristles and add an approximately pea-sized amount of toothpaste to the bristle-head
- Distribute toothpaste evenly onto your teeth, while no brushing necessary.
- Turn device on, 1 to 3 green LEDs will turn on sequentially and the device will vibrate gently indicating that the Ultrasound wave cleaning process started.
- Place tooth bristle-head gently against teeth and gums. Ultrasound waves and vibrations will gently. Do NOT MOVE the tooth bristle-head, no brushing required. Do NOT PRESS against teeth or gums.
- After approximately 6 seconds the vibrations will shortly stop to indicate you should move to the next cleaning section.
- Treat all teeth and gum surface areas sequentially with Ultrasound, no brushing necessary
- When all teeth/gums surfaces have been cleaned – gently move bristles over teeth and gums to remove loosened up bacterial plaque.
- With device still turned on - rinse bristle-head under running water – this will clean away bacteria from the bristles.
- Rinse your mouth – and you are done.



- Batteries last up to 2 weeks** with 2x3 minutes daily usage. Device needs to be charged when only 1 green LED lights up. For best battery maintenance keep the device in the charging unit.
- For best cleaning results - replace bristle-head at least every 2-3 months.** Ultrasound waves will weaken bristles and over time resulting in a lower cleaning efficiency.

For more information

- Check our website – www.emmident-ultrasound.com
- Check our emmi-dent video instructions at <https://youtu.be/heb3-8Qlqq5>
- Contact us at CustomerService@emmident-ultrasound.com

Distributed by Ultra Oral Care, Inc., USA

