



## Ultrasound Toothpastes with Nano-Bubbles Technology

These specially formulated Emmi-dent Ultrasound Toothpastes, used together with Emmi-dent's original Ultrasound device, will ensure best cleaning results. The formula is enriched with millions nano-cleaning bubbles that penetrate even the smallest dental crevices and gaps, removing impurities, traces of food, germs and harmful bacterial plaque (= "cavitation" cleaning process).

The toothpaste is part of emmi-dent's 100% Ultrasound technology cleaning system and cleans teeth and gums extra gentle and without any abrasion because it cleans without brushing. Emmi-dent toothpastes do NOT contain abrasive micro-particles. Using Emmi-dent Ultrasound your teeth will not lose their natural protection against bacteria and the risk of dental diseases decreases.

Even people with most sensitive, inflamed or receding gums can confidently and pain-free clean their teeth and gums, knowing that emmi-dent Ultrasound will deep-clean and gently remove harmful bacterial plaque. Emmi-dent Ultrasound toothpastes are available in the following variants:

- **Fresh** - fresh mint taste with fluoride.
- **Whitening** - extra fresh mint taste and extra fluoride.
- **Nature with Calcium** - mild taste, non-mint flavor, does not contain fluoride.
- **Nature** - extra mild/neutral taste, non-mint flavor, does not contain fluoride.
- **Children** - fruity flavor, with reduced fluoride content.

**DIRECTIONS** – adults and children 12 yrs. of age and older: Clean teeth thoroughly, preferably after each meal or at least twice a day, or as directed by a dentist or physician. Children under 12 yrs. of age: ask a dentist.

**WARNING** – Keep out of reach for children under 12 yrs. of age. If more than used for brushing is accidentally swallowed, get medical help or contact a Poison Control Center right away.

**SHELF LIFE** - Minimum shelf life is 3 years from the manufacturing date. Depending on the printed code format you can read the manufacturing year as follows: Format "JTXXY" - the first digit "J" stands for the year (e.g. 9 = 2019) or Format "MM.JJ/XY" – the 3rd and 4th digits "JJ" stand for the year (e.g. 19 = 2019)

### **emmi®-dent Toothpastes' Ingredients**

**Fresh:** • Water • Hydrated Silica • Sorbitol • Propylene Glycol • Tetrapotassium Pyrophosphate • Xanthan Gum • Sodium C14-16 Olefin Sulfonate • Mentha Arvensis Leaf Oil – Aroma • Sodium Fluoride • CI 77891 (Titanium Dioxide) • Sodium Saccharin • Allantoin • Chamomilla Recutita Flower Extract – Aroma • Salvia Triloba Leaf Extract – Aroma • Sodium Methylparaben. Fresh does NOT contain gluten, nanomaterials or micro-plastics.

AROMA-percentage: > 1% < 5% / FLUORIDE: < 0.1% / MENTHA ARVENSIS LEAF OIL: > 1% < 5%

**Nature:** • Aqua • Hydrated Silica • Sorbitol • Propylene Glycol – Propylenglykol • Sodium C14-16 Olefin Sulfonate – Natrium C14-16 Olefin Sulfonat • Aroma • CI 77891 (Titanium Dioxide, E171) • Cellulose Gum – Cellulose • Stevia Rebaudiana Extract. Does NOT contain fluoride, gluten, nano-materials or micro-plastics. This product is vegan.

AROMA-percentage: > 1% - < 5%

**Nature with calcium:** • Water • Calcium Glycerophosphate • Hydrated Silica • Sorbitol • Propylene Glycol • Sodium C14-16 Olefin Sulfonate • Aroma • CI 77891 (Titanium Dioxide) • Cellulose Gum • Stevia Rebaudiana Extract – Aroma Nature Calcium is vegan, does NOT contain fluoride, gluten, nano materials or micro plastics. AROMA: > 1% < 5%

**Whitening:** • Aqua • Sorbitol • Hydrated Silica • Sodium Bicarbonat – Natrium Hydrogencarbonat • Propylene Glycol – Propylen-Glykol • Pentasodium Triphosphate – Pentanatrium-Triphosphat • Sodium C14-16 Olefin Sulfonate – Natrium C14-16 Olefin Sulfonat • Cellulose Gum (E 466) • Aroma - Aroma • Sodium Fluoride • Zinc Chloride • Sodium Saccharin – Natrium-Saccharin • Limonene • CI 74160 (Copper Phthalocyanine) - • CI 77891 (Titanium Dioxide, E171). Whitening does NOT contain gluten, nanomaterials or micro-plastics. AROMA-percentage: > 1% < 5% / FLUORIDE: < 0.1%

**Kids (with reduced fluoride content):** • Water • Hydrated Silica • Sorbitol • Propylene Glycol • Sodium C14-16 Olefin Sulfonate • Cellulose Gum • Sodium Fluoride\* • Sodium Saccharin • Sodium Methylparaben • Titanium Dioxide (CI 77891) • Aroma Kids does NOT contain gluten, nano materials or micro plastics. \*) contains 500 ppm (parts per million) fluoride.



## Bristle-Head Attachments with patented, embedded Piezo-Chip Technology

The emmi-dent Oral Waver bristle-head attachment series is part of its patented 100% Ultrasound oral health cleaning system.

These bristle-head attachments work with 100% Ultrasound wave technology and gently deep-clean without any abrasive brushing.

The ultrasound waves are generated by the patented Piezo-Chip Technology that is embedded inside the bristle-head. The Piezo-Chip generates ultrasound waves with up to 96 Million oscillations per minute. The specially designed bristles then serve as antennas emitting the ultrasound waves for most effective deep-cleaning of teeth and gums.

The combination of the emmi-dent ultrasound waves and the ultrasound nano-bubble toothpaste gently deep-clean teeth and gums even inside tiniest crevices where ordinary toothbrushes cannot reach – all without abrasive brushing.

Emmi-dent's 100% ultrasound wave technology – is suited for all and especially recommended for people with sensitive teeth and gums, periodontitis, inflamed gums, braces, implants and after oral surgery.

Emmi-dent offers a range of different bristle-head designs for different oral health applications:

1. **Regular** bristle-head size and bristle design
2. **Wave** bristle-head designed for better inter-dental cleaning
3. **Braces** - bristle-head designed for best cleaning of braces
4. **Kids/Young adults** – with a smaller sized bristle-head

### Check for the correct bristle-head attachment when ordering!

Bristle-head attachments are available in 2 different connection designs:

- **4-PIN connection for Platinum and Oral Waver devices**
- **2-PIN connection for Chrome/Metallic and Professional devices**

When ordering please check whether you need 2-PIN or 4-PIN attachments to make sure that the attachments that can be used with your emmi-dent device!

Questions? - Contact us at: [CustomerService@emmident-ultrasound.com](mailto:CustomerService@emmident-ultrasound.com)  
To reorder and special offers check our website: <https://www.emmident-ultrasound.com>